


## Direct assessment of the availability of iron for haemoglobinisation

Functional iron deficiency (FID) is a state in which iron incorporation into erythroid precursor cells is insufficient despite apparently normal body iron stores. The most common test to determine functional iron availability is the assessment of transferrin saturation (TSAT). However, TSAT may be affected by factors that are unrelated to iron status, such as infection and inflammation. A sensitive and easily accessible blood test marker would therefore be desirable to rapidly evaluate the iron that is directly available for haemoglobin synthesis. There is a parameter backed with evidence from studies that can be used: the reticulocyte haemoglobin equivalent (RET-He).



A 68-year-old female with symptoms of fatigue and difficulty in coping with activities of daily life presented to her nephrologist. She was suffering from chronic kidney disease at stage III due to polycystic kidney disease. The blood count revealed anaemia with a decreased haemoglobin (9.9 g/dL or 6.14 mmol/L) and reticulocyte haemoglobin equivalent (27.4 pg or 1.70 fmol). Due to the underlying inflammation, traditional tests such as serum iron (36 µg/dL), ferritin (788 ng/mL) and TSAT (15%) were not fully conclusive but indicated functional iron deficiency. The patient was referred for the administration of erythropoiesis-stimulating agents (ESA) to promote the production of red blood cells. Supplemental iron was given prior to ESA therapy to ensure adequate haemoglobinisation of the new red blood cells. The patient completed five weeks of I.V. iron therapy with good response already one week after the therapy was initiated, with RET-He levels increasing to 29.8 pg or 1.85 fmol. This increase in RET-He indicated the improved availability of iron for haemoglobinisation and corroborated the effectiveness of treatment [1]. This was further substantiated by the patient's reduction in symptoms.

### What is the reticulocyte haemoglobin equivalent, or RET-He?


- RET-He is a haematology parameter which reflects the haemoglobin content of reticulocytes – immature red blood cells.
- RET-He provides an early assessment of the available iron that was utilised in the red blood cell production over the previous 2–4 days [2, 3].
- RET-He reference range: 29.3–35.4 pg or 1.82–2.20 fmol [4].
- RET-He has been reported to have high accuracy, sensitivity and specificity for identifying iron deficiency [5, 6].
- The test methodology is based on fluorescence flow cytometry.
- RET-He is readily available from a routine EDTA blood sample analysis in the laboratory.

RETICULOCYTE  
HAEMOGLOBIN  
CONTENT  
CLINICAL USE

Know more.

Decide with confidence.

Act faster.



## Benefits of RET-He with functional iron deficiency

- Assesses the content of haemoglobin in reticulocytes
- Reflects the bioavailability of iron for erythropoiesis, comparable to transferrin saturation [2, 11]
- Not affected by the acute-phase reaction [7, 8]
- Affected by biological variation to a much lower degree than TSAT and ferritin [9]
- Values below 29 pg or 1.80 fmol are indicative of functional iron deficiency in the context of postoperative anaemia and anaemia in chronic kidney disease (CKD) [10, 11]
- Early indicator of the response to iron therapy and/or ESA within 2–4 days [1, 12]

Benefit from more background information in our freely accessible white papers:

[www.sysmex-europe.com/whitepapers](http://www.sysmex-europe.com/whitepapers)

Accessing our literature lists will give you an overview of more peer reviewed articles on this topic:

[www.sysmex-europe.com/publications](http://www.sysmex-europe.com/publications)

- [1] **Almashjary MN et al. (2022):** Reticulocyte Hemoglobin-Equivalent Potentially Detects, Diagnoses and Discriminates between Stages of Iron Deficiency with High Sensitivity and Specificity. *J Clin Med.* 2022 Sep 26; 11(19): 5675.
- [2] **Hoennemann C et al. (2021):** Reticulocyte and Erythrocyte Hemoglobin Parameters for Iron Deficiency and Anemia Diagnostics in Patient Blood Management. A Narrative Review. *J Clin Med.* 2021 Sep; 10(18): 4250.
- [3] **Hoennemann C et al. (2021):** Reticulocyte Haemoglobin as a Routine Parameter in Preoperative Iron Deficiency Assessment. *Endocrinol Metab Vol. 5 No.1:* 154.
- [4] **van Pelt JL et al. (2022):** Reference intervals for Sysmex XN hematological parameters as assessed in the Dutch Lifelines cohort. *Clin Chem Lab Med.* 2022 May 2; 60(6): 907–920.
- [5] **Buttarelo M et al. (2016):** Evaluation of the hypochromic erythrocyte and reticulocyte hemoglobin content provided by the Sysmex XE-5000 analyzer in diagnosis of iron deficiency erythropoiesis. *Clin Chem Lab Med.* 2016 Dec 1; 54(12): 1939–1945.
- [6] **Mehta S et al. (2016):** Reticulocyte Hemoglobin vis-a-vis Serum Ferritin as a Marker of Bone Marrow Iron Store in Iron Deficiency Anemia. *J Assoc Physicians India.* 2016 Nov; 64(11): 38–42.
- [7] **Thomas L et al. (2005):** Reticulocyte hemoglobin measurement-comparison of two methods in the diagnosis of iron-restricted erythropoiesis. *Clin Chem Lab Med.* 2005; 43(11): 1193–202.
- [8] **Thomas C et al. (2002):** Biochemical markers and hematologic indices in the diagnosis of functional iron deficiency. *Clin Chem.* 2002 Jul; 48(7): 1066–76.
- [9] **Van Wyck DB et al. (2010):** Analytical and biological variation in measures of anemia and iron status in patients treated with maintenance hemodialysis. *Am J Kidney Dis.* 2010 Sep; 56(3): 540–6.
- [10] **Shander A et al. (2022):** Recommendations from the International Consensus Conference on Anemia Management in Surgical Patients (ICCAMS). *Ann Surg.* 2022 Sep 21; Online ahead of print.
- [11] **Muñoz M et al. (2017):** International consensus statement on the peri-operative management of anaemia and iron deficiency. *Anaesthesia.* 2017 Feb; 72(2): 233–247.
- [12] **Auerbach M et al. (2021):** Using Reticulocyte Hemoglobin Equivalent as a Marker for Iron Deficiency and Responsiveness to Iron Therapy. *Mayo Clin Proc.* 2021 Jun; 96(6): 1510–1519.